When the cover of the device is closed (in vaccination mode) pay attention to the needle length to be 18-22 mm depending on chick size.

The space between the vaccination button and needle should be approximately 1-2 mm. Otherwise there may be bleeding or out dispensing during the vaccination.

The space between the needle tip and treatment block's upper side should be 30-34 mm.

Check the position of the needle according to these instructions before each vaccination run.

It is recommended to replace the needle after every 2000 vaccination.

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BEFORE BEGINNING THE VACCINATION CHECK THESE:

- ✓ The operator should follow safety guidelines,
- ✓ The vaccine should be connected to the injectors and dosage should be adjusted,
- ✓ Compressed air pressure should be 4-6 bars,
- ✓ Injectors should be checked using the positioning switches,
- ✓ Needle's position and direction should be as instructed,
- ✓ Counter should be set to target number.

If above points are met then the device is ready for operation.

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4. Cleaning and Maintenance

4.1. Air Supply and Filter


- Make sure the air supplied is dry and clean.
- Pay attention to the air filter and dust container to be clean.

To clean the dust container and filter:

- First disconnect air supply, then turn the filter to the left. Remove and clean the filter. Replace the filter by turning to the right.
- Restart the device and check.

Just press the part #4 under the container to drain the water.

4.2. Fixed and Adjustable Injectors

- The injector compartment under the cover is simply composed of two parts, fixed and moving.
- Perform daily cleaning using compressed air and alcohol.
- It is recommended to use only dampened cloth to wipe and not use water jet or water with detergent.

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
4.3. Treatment Block

- This part should be cleaned after each vaccine change or daily with compressed air and alcohol.


- For weekly cleaning position adjustment screw should be unscrewed, the bottom of the treatment block and the area under it should be wiped using a dampened cloth, and then cleaned with compressed air and alcohol.

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4.4. Injector Plate

- The plate should be wiped using a dampened cloth, and then cleaned with compressed air and alcohol, vaccine residues should be removed at the end of every working day.

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4.5. Injectors

- No vaccine should remain in the injectors and check valves after vaccination is completed, so clean the injectors by drawing in first alcohol and then water.
- Make sure there is no vaccine or alcohol left in the injectors.
- After you cleaned the check valves with alcohol and water you should separate and clean vaccine container and injector axis.
- You can use alcohol and water for this, alternatively you can keep the parts in 60°C water for 10 minutes.
- Make sure there is no vaccine and alcohol left in the injectors.
- You can dry the injectors using a clean cloth or compressed air.

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PERFORM DAILY

- Clean the inner and outer parts of the device.
- Keep the treatment block clean.
- Clean the vaccination button on the treatment block.
- Keep moving and fixed parts clean.
- Make sure the injectors are clean and not contaminated.
- Perform cleaning activities using compressed air and dampened cloth.

PERFORM PERIODICALLY (DEPENDING ON THE FREQUENCY OF USAGE)

- Make sure supplied air and air filter is clean.
- Check there is sufficient amount of oil.
- Clean the bottom of the treatment block.
- Clean the injector plate.
- Ask technical personnel to perform monthly maintenance.

OPERATE THE DEVICE IN A DRY ENVIRONMENT,
STORE AT ROOM TEMPERATURE, PROTECT FROM
DIRECT SUNLIGHT

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5. Periodic Replacement of Parts

- O-rings: Every 250,000 vaccinations
- Needle: Every 2,000 vaccinations
- Check valve springs (#9, #10): Every 250,000 vaccinations
- In case of faulty parts contact technical service
- Replace parts only with original ones supplied by Novakim Mectech.

6. Technical Points

- Vaccination options:
 - Water based (with adjustable injector)
 - Water based (with fixed injector) and Oil based (with adjustable injector)
 - Oil based (with fixed injector) and water based (with adjustable injector)
 - Oil based (with adjustable injector)

You can assign priority to a vaccine by connecting the vaccine you want apply first to the adjustable injector, and the second one to the fixed injector.
- 3000-4500 vaccinations per hour can be performed regardless of single or double vaccination.
- The device counts the number of vaccination in two ways; daily counter and application counter.
- For 2000 doses of vaccination, deviations are ± 5 doses for oil based vaccines and ± 14 doses for water based vaccines.

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7. Possible Problems and Alternative Solutions

- In this section alternative solutions to possible problems which may arise during installation or operation of the device are given:
- Steps to follow in order are described under each problem.
- For each problem follow in order the steps described.
- Do not open the service cover of the device and intervene the pneumatic parts.
- Contact technical service if the problem is not solved.
- Use only original parts supplied by Novakim Mectech to replace old or faulty parts.
- Do not neglect daily and periodic cleaning and maintenance.

IN THE BEGINNING

➤ **Pressure gauge shows 0 bar:**

- Check the air pressure from the compressor.
- Check tubing connection.
- Make sure there is no air leakage.
- If you detect air leakage find its source and stop leakage.

➤ **Pressure gauge shows below 4 bars:**

- Check the air pressure from the compressor.
- Check tubing connection.
- Make sure there is no air leakage.
- If you detect air leakage find its source and stop leakage.
- Adjust pressure to desired value.
- If the pressure still won't go up ask technical personnel to check the compressor.

➤ **Pressure gauge shows above 6 bars:**

- Try to lower the pressure with pressure adjustment knob.
- Contact technical service if pressure can't be lowered.

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- Check connections in the treatment block.
- If you detect air leakage find its source and stop leakage.
- **Vaccination button won't start the operation:**
 - Press RESET button.
 - Clean treatment block.
 - Make sure the counter is set to value other than 000.
 - If the vaccination button won't move take it out, clean and place it back.

DURING VACCINATION

- **Chick has vaccine residue on it:**
 - Check the way the operator works.
 - Inspect check valve system of the injectors and replace necessary parts.
 - Check needle depth (18-22 mm)
 - Check needle's position and direction.
 - Check the space between the needle and vaccination button (1-2 mm).
 - Check the position of the treatment block.
- **Chicks bleed:**
 - Check the way the operator works.
 - Check needle depth (18-22 mm)
 - Check needle's position and direction.
 - Check needle's tip.
 - Check the space between the needle and vaccination button (1-2 mm).
 - Check the position of the treatment block.
- **Vaccination button won't start operation:**
 - Press RESET button.
 - Clean the treatment block.
 - Make sure the counter is set to value other than 000.
 - If the vaccination button won't move take it out, clean and place it back.

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8. List of Device's Components and Spare Parts

PVT VACCINATOR

PVT-K
PVT-HB
PVT-SB
PVT-ES
PVT-CH

PVT-CH

PVT-CH-01
PVT-CH-02
PVT-CH-03
PVT-CH-04
PVT-CH-05
PVT-CH-06
PVT-CH-07
PVT-CH-08
PVT-CH-09
PVT-CH-10
PVT-CH-11
PVT-CH-12
PVT-CH-13
PVT-CH-14
PVT-CH-15
PVT-CH-16
PVT-CH-17
PVT-CH-18
PVT-CH-19
PVT-CH-20

PVT-ES

PVT-ES-01
PVT-ES-02
PVT-ES-03
PVT-ES-04
PVT-ES-05

PVT-HB

PVT-HB-01
PVT-HB-02
PVT-HB-03
PVT-HB-04
PVT-HB-05
PVT-HB-06
PVT-HB-07
PVT-HB-08
PVT-HB-09
PVT-HB-10
PVT-HB-11
PVT-HB-12
PVT-HB-13
PVT-HB-14
PVT-HB-15
PVT-HB-16
PVT-HB-17
PVT-HB-18
PVT-E-1

PVT-SB

PVT-SB-01
PVT-SB-02
PVT-SB-03
PVT-SB-04
PVT-SB-05
PVT-SB-06
PVT-SB-07
PVT-SB-08
PVT-SB-09
PVT-SB-10

COVER

PVT-K-01
PVT-K-02
PVT-K-03
PVT-K-04

PVT-K-05
PVT-K-06
PVT-K-07
PVT-K-08
PVT-K-09
PVT-K-10
PVT-K-11
PVT-K-12

OIL BASED INJECTOR

PVT-E-001
PVT-E-002
PVT-E-003
PVT-E-004
PVT-E-005
PVT-E-006
PVT-E-008
PVT-E-010
PVT-E-011
PVT-E-012
PVT-E-013
PVT-E-016

WATER BASED INJECTOR

PVT-E-001
PVT-E-002
PVT-E-003
PVT-E-004
PVT-E-005
PVT-E-006
PVT-E-007
PVT-E-008
PVT-E-009
PVT-E-010
PVT-E-011
PVT-E-012
PVT-E-015

PVT VACCINATOR

PVT-K
PVT-HB
PVT-SB
PVT-ES
PVT-CH

PVT-CH

PVT-CH-01
PVT-CH-02
PVT-CH-03
PVT-CH-04
PVT-CH-05
PVT-CH-06
PVT-CH-07
PVT-CH-08
PVT-CH-09
PVT-CH-10
PVT-CH-11
PVT-CH-12
PVT-CH-13
PVT-CH-14
PVT-CH-15
PVT-CH-16
PVT-CH-17
PVT-CH-18
PVT-CH-19
PVT-CH-20

PVT-ES

PVT-ES-01
PVT-ES-02
PVT-ES-03
PVT-ES-04
PVT-ES-05

PVT-HB

PVT-HB-01
PVT-HB-02
PVT-HB-03
PVT-HB-04
PVT-HB-05
PVT-HB-06
PVT-HB-07
PVT-HB-08
PVT-HB-09
PVT-HB-10
PVT-HB-11
PVT-HB-12
PVT-HB-13
PVT-HB-14
PVT-HB-15
PVT-HB-16
PVT-HB-17
PVT-HB-18
PVT-E-1

PVT-SB

PVT-SB-01
PVT-SB-02
PVT-SB-03
PVT-SB-04
PVT-SB-05
PVT-SB-06
PVT-SB-07
PVT-SB-08
PVT-SB-09
PVT-SB-10

COVER

PVT-K-01
PVT-K-02
PVT-K-03
PVT-K-04
PVT-K-05
PVT-K-06
PVT-K-07
PVT-K-08
PVT-K-09
PVT-K-10
PVT-K-11
PVT-K-12


OIL BASED INJECTOR

PVT-E-001
PVT-E-002
PVT-E-003
PVT-E-004
PVT-E-005
PVT-E-006
PVT-E-008
PVT-E-010
PVT-E-011
PVT-E-012
PVT-E-013
PVT-E-016

**WATER BASED
INJECTOR**

PVT-E-001
PVT-E-002
PVT-E-003
PVT-E-004
PVT-E-005
PVT-E-006
PVT-E-007
PVT-E-008
PVT-E-009
PVT-E-010
PVT-E-011
PVT-E-012
PVT-E-015

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



Unit 6, No. 101, West Bagherkhan, Befor Tohid Sq.
Chamran Highway.
Tehran, Iran
Postal code: 1441415617
Tet: +98(21)66935254 ,
Fax: +98(21) 66592064
Parsian Pakhsh exir


MOVING PART


OILY INJECTOR

WATER INJECTOR

PVT-SE

COVER